Pulido
ELECTROLITICO

R-080490

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 Practical, efficient, designed for use in dental mechanics. Up to three frameworks can be polished at once.

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
OPERATING PRINCIPLE

Under a microscope, all metals show a more or less rough surface with protrusions of all types. If they are placed in an electrolytic bath (fig. 1) and connected to an anode a curious thing happens: the electrolyte dissociates itself, producing charged molecules called ions which literally pull atoms off the surface of the metal and

combine with them. The particles produced in this way are electrically charged, and the difference in potential between the cathode and anode (fig. 2) moves them gradually towards the cathode. The particles of metal are deposited on the cathode (fig. 3), forming more charged ions, and the cycle starts over. The process ends only when the electricity is cut off.

The areas on the metal with the biggest protrusions have more surface to come into contact with the electrolyte and therefore lose most atoms during polishing. As a result when polishing is completed the surface of the metal is much less rough (fig. 4), and the protrusions are much smoother. This enables it to reflect light more evenly, giving the metal more shine. The excess metal removed from the anode (+) ends up deposited on the surface of the cathode.

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
Rounded frame made entirely of stainless steel.

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Tank and lid made of ABS injected thermoplastic. Being an independent part, it is easy and safer to handle for cleaning and maintenance.

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