Instructions for Use

TF™ Adaptive

Rotary and Hand Files

Kerr Italia S.r.l.
CE 2797 Via Passanti, 174
Scafati (SA), 84018 Italy
+39 081 850 8311
kerrdental.com

UK CA 0086 Authorised Rep in UK:
Kerr UK Ltd.
Orega Stockley Park
4 Longwalk Road
Stockley Park
Uxbridge UB 11 1FE
United Kingdom

077-0728e_Rev.00
2026-06

Rx only* *for US law only | NON STERILE | MD | 132°C 270°F | 121°C 250°F

Kerr™ ENDODONTICS

©2026 Kerr Corporation. All Rights Reserved.

SYMBOLS DESCRIPTION

Full explanation of symbols used on Kerr packaging is located at:
http://www.kerrdental.com/symbols-glossary

Rx only* *for US law only	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.
121°C 250°F	Gravity Displacement Sterilizable in a steam sterilizer (autoclave) at 121°C (250°F).
132°C 270°F	Pre-Vacuum Sterilizable in a steam sterilizer (autoclave) at 132°C (270°F).

en - ENGLISH

TF™ ADAPTIVE

Rotary and Hand Files

GENERAL INFORMATION

Rotary files and hand files are used for root canal shaping and cleaning during endodontic treatment. Files are offered in either refill or assortment packs.

Intended Purpose: Rotary and hand files are used for root canal shaping and cleaning during endodontic treatment.

Intended Users: This product is to be used by qualified dental professionals.

Intended Patient Population: The device is for use in a general population of patients in need of a root canal treatment.

Clinical Benefits: It helps to clean the root canal walls so the endodontically treated tooth can be restored and remain in function.

PRODUCT COMPOSITION

Rotary files contain Brass, Stainless Steel, Nickel-Titanium, and a latex-free silicone stopper. Hand files contain Stainless Steel, Polyetherimide, and a latex-free silicone stopper.

INDICATIONS FOR USE

Rotary and hand files are instruments used for shaping, cleaning, and enlarging root canals for injured pulp conditions.

CONTRAINDICATIONS

None known.

WARNINGS

- Rotary files are single use products and must be discarded after each use. Do not re-use. Risk of cross-contamination and breakage.
- All files are non-sterile and must be cleaned and sterilized prior to use.
- These products should not be used for individuals with sensitivity to brass, nickel-titanium or stainless steel. The dental professional should inform the patient that this product is known to cause a brass, nickel-titanium and stainless steel allergic reaction.
- Use of a rubber dam is recommended when using these instruments for root canal procedures to avoid accidental aspiration or ingestion.
- In canals with severe or sudden apical curvature additional caution should be exercised due to the heightened risk of file separation.

PRECAUTIONS

- Determining the reference working length with the Endo Stop is imperative to ensure proper instrumentation.
- Rotary files are made to be used with a slow speed handpiece. Ensure that the rotary file is securely attached to the handpiece prior to use.
- When using the hand files ensure straight-line access and do not force a file down the canal as it might separate, get stuck or create transportation.
- TF Adaptive files are designed to be used with elements™ e-motion motor.

ADVERSE EVENTS

If a serious incident occurs with this medical device, report it to the legal manufacturer, Kerr Italia, S.r.l., and to the competent authority for that country in which the user and/or patient is established.

STEP-BY-STEP INSTRUCTIONS

Limitations of Use:

- Hand files can be used a maximum of five (5) times.
- Rotary files are single-use only and should be discarded after each use.
- Extremely curved canals. In case they need to be used, use a new file per canal.
- When the file shows the following:
 - Untwisted threads
 - Dull cutting surfaces
 - Bent instrument
 - Signs of corrosion
 - Unable to sterilize the instrument prior to use

TF Adaptive Instrumentation Technique:

1. Locate the canal orifices using straight line access and obtain patency with the hand K-Files.
2. Begin the instrumentation technique by taking a coronal shaper to resistance.
3. Establish working length using K-Files, apex locator, and radiograph as needed.
4. Use the "TF Adaptive" setting on your elements e-motion motor.
5. Canal size and file sequence determination (Refer to Table 1. for File Size):
 - Small Canals (SM) – Using tactile feel, if you cannot advance using a #15 K-File to working length, then the canal is Small. Use the Small Pack (one color band) and instrument in sequence as indicated below.
 - Medium / Large Canals (ML) - Using tactile feel, if a #15 K-File is loose, then the canal size is Medium/Large. Use the Medium/Large Pack (two color bands) and instrument in sequence as indicated below.
6. Ensure the pulp chamber is filled with NaOCl (Sodium Hypochlorite), or 17% EDTA and make sure the file is rotating as you engage the canal.
7. Slowly advance the green (SM1 or ML1) with a single controlled motion until the file engages dentin and then completely withdraw the file from the canal. Do not force apically. Do not peck.
8. Wipe off the flutes. Deliver irrigant to the pulp chamber and confirm canal patency with a #15 K-File.
9. Repeat steps 7 and 8 using the file you started with until working length is achieved.
10. Repeat steps 7 and 8 with the yellow SM2 or ML2 until the file reaches working length. If the desired apical size is achieved the sequence is complete.
11. For larger apical sizes, repeat steps 7 and 8 with the red SM3 or ML3 until the file reaches working length.
12. Obturate: Follow preferred Obturation technique.

Table 1. - File Size Reference Chart:

SMALL (SM)	SM1: #20/ .04	SM2: #25/ .06	SM3: #35/ .04
MEDIUM/ LARGE (ML)	ML1: #25/ .08	ML2: #35/ .06	ML3: #50/ .04

Note:

Techniques suggested are guidelines for the average case. It is imperative that clinical judgment always be exercised, appropriate setting adjustments are made, and additional safety measures are implemented based on individual case requirements.

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