Instructions for Use

ZenFlex™

Kerr™ ENDODONTICS

Kerr Italia S.r.l.
Via Passanti, 174
2797 Scafati (SA), 84018 Italy
+39 081 850 8311
kerrdental.com

Distribution Facility:
Kerr Australia Pty. Ltd.
Unit 6, 12 Mars Road
Lane Cove West, NSW 2066
Australia
+61 2 8870 3000

UK CA 0086 Authorised Rep in UK:
Kerr UK Ltd.
Orega Stockley Park
4 Longwalk Road
Stockley Park
Uxbridge UB 11 1FE
United Kingdom

077-0648e_Rev.00
2026-06

©2026 Kerr Corporation. All Rights Reserved.

SYMBOLS DESCRIPTION

Full explanation of symbols used on Kerr packaging is located at:
http://www.kerrdental.com/symbols-glossary

	Gravity Displacement Sterilizable in a steam sterilizer (autoclave) at 121°C (250°F).
	Pre-Vacuum Sterilizable in a steam sterilizer (autoclave) at 132°C (270°F).

en - ENGLISH

ZenFlex™
Rotary Files

GENERAL INFORMATION

ZenFlex rotary files are used for root canal shaping and cleaning during endodontic treatment in rotary motion.

Files are offered in either 0.4 and 0.6 tapers and as a refill or assortment packs of 20, 25, 30, 35, 40, and 45 tips. Additional refills are available in 50 and 55 tips. ZenFlex files are available in 21, 25, and 31 mm in length.

Intended Purpose: ZenFlex files are used for root canal cleaning and shaping during endodontic treatment.
Intended Users: This product is to be used by qualified dental professionals.
Intended Patient Population: The device is for use in a general patient population in need of a root canal treatment.
Clinical Benefit: Allows the dentist to efficiently prepare and shape the root canal for cleaning and obturation during endodontic treatment.

PRODUCT COMPOSITION

These products contain a brass handle, nickel-titanium alloy, and a latex-free silicone stopper.

INDICATIONS FOR USE

ZenFlex rotary files are instruments used for shaping, cleaning, and enlarging root canals for injured pulp conditions.

CONTRAINDICATIONS

None known.

WARNINGS

- These products should not be used for individuals with sensitivity to nickel or titanium. Dentists shall inform the patient that nickel and titanium are known to cause an allergic reaction.
- ZenFlex is a single-use product. Discard after use. Do not reuse as there is a risk of cross-contamination and breakage.
- ZenFlex is a non-sterile, single-use product and must be cleaned and sterilized prior to use.
- Use of a rubber dam is recommended when using these instruments for root canal procedures to avoid accidental aspiration or ingestion.
- In canals with severe or sudden apical curvature additional caution should be exercised due to the heightened risk of file separation.
- ZenFlex files should not be used at a speed and/or torque higher than what is recommended by the manufacturer as the risk of apical perforation, apical transportation, strip perforation, file separation, and/or root fracture is possible.

PRECAUTIONS

- Determining the reference working length with the Endo Stop is imperative to ensure proper instrumentation.
- ZenFlex is made to be used with an Endodontic handpiece. Ensure that the rotary file is securely attached to the handpiece prior to use. Refer to Table 1 for torque and speed recommendations.
- Use personal protective equipment such as gloves and safety glasses.
- Ensure that the rotary file is securely attached to the handpiece prior to use. Pull the file to make sure that the file is properly latched in the contra-angle.

ADVERSE EVENTS

If a serious incident occurs with this medical device, report it to the legal manufacturer, Kerr Italia, S.r.l., and to the competent authority for that country in which the user and/or patient is established.

STEP-BY-STEP INSTRUCTIONS

Limitations on Use:

- ZenFlex files are single-use only and should be discarded after each use.
- Extremely curved canals. In case they need to be used, use a new file per canal.
- Unable to sterilize the instrument prior to use.
- When the file shows the following:
 - Untwisted threads
 - Dull cutting surfaces
 - Bent instrument
 - Signs of corrosion

Rotary Files Technique:

1. Locate the canal orifices using straight-line access and obtain patency with hand files.
2. Establish working length using hand files, apex locator, and radiograph as needed.
3. Establish the glide path using hand files or rotary glide path files to minimize file separation, ledging and transportation.
4. Take the first rotary file in the sequence (20 tip) to resistance, repeating as needed for difficult, curved, or narrow canals.
5. Take each instrument to resistance at the torque and speed recommended in the Table 1. for no longer than 5-7 seconds.
6. Additional Clinical Tips:
 - Irrigate the canal using NaOCl (Sodium Hypochlorite), 17% EDTA or a suitable lubricant agent to facilitate debris removal.
 - Use light pressure and avoid force. When the instrument no longer advances apically, proceed to a smaller file.
 - Prior to use, check instruments for defects; including deformation, breakage, fractures, corrosion, handle damage, and loss of color coding or size marking. Discard any product with defects.

Shaping Technique

1. Use settings for your handpiece according to Table 1. below.
2. With the rotary file rotating as you enter the canal, slowly advance the rotary file with a single controlled motion until the file engages dentin, then completely withdraw the rotary file from the canal. Do not force the rotary file apically. Do not peck.
3. Wipe the flutes. Irrigate with NaOCl then confirm patency with #10 K-File.
4. Repeat step 2 using the rotary file you started with until you are not able to passively advance or working length is achieved. Never apply pressure on the rotary file.
5. Repeat with the next rotary file in chosen sequence until the final working length and desired shape are achieved.

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