Instructions for Use

K-Files™, K-Flex™, Triple-Flex™, Hedstrom™, and Reamers


Hand Files

Kerr Italia S.r.l.
CE 2797 Via Passanti, 174
Scafati (SA), 84018 Italy
+39 081 850 8311
kerrdental.com

Distribution Facility:
Kerr Australia Pty. Ltd.
Unit 6, 12 Mars Road
Lane Cove West
NSW 2066, Australia
+61 2 8870 3000

077-1129e_Rev.00
2026-06

UK CA 0086 Authorised Rep in UK:
Kerr UK Ltd.
Orega Stockley Park
4 Longwalk Road
Stockley Park
Uxbridge UB 11 1FE
United Kingdom

Kerr™ ENDODONTICS

©2026 Kerr Corporation. All Rights Reserved.

SYMBOLS DESCRIPTION

Full explanation of symbols used on Kerr packaging is located at:
http://www.kerrdental.com/symbols-glossary

121°C (250°F)	Gravity displacement Sterilizable in a steam sterilizer (autoclave) at 121 °C (250 °F).
132°C (270°F)	Pre-Vacuum Sterilizable in a steam sterilizer (autoclave) at 132 °C (270 °F).
134°C (273°F)	Pre-Vacuum Sterilizable in a steam sterilizer (autoclave) at 134 °C (273 °F).

en - ENGLISH

K-Files™, K-Flex™, Triple-Flex™, Hedstrom™, and Reamers
Hand Files

GENERAL INFORMATION

Endodontic hand files are used for root canal cleaning and shaping during endodontic treatment.

K-Files, K-Flex, and **Triple-Flex** files are used for the negotiation of canals and for the initial enlargement of canals to make way for rotary instrumentation increasing flexibility to negotiate with curved canals and stressful situations.
Hedstrom files are used for rapid enlargement of canals because of the aggressive cutting. They are also used for the removal of root canal filling materials e.g., gutta percha during retreatment of the root canal.
Reamers are used for widening the root canals to clean and shape canals using a rotary motion providing flexibility with integrity.

All hand files are available in lengths of 21, 25 and 30 mm, in packets of 6 files and in sizes ranging from 06 to 80.

Intended Purpose: Endodontic hand files are used for root canal cleaning and shaping during endodontic treatment.
Intended Users: This product is to be used by qualified dental professionals.
Intended Patient Population: The device is for use in a general population of patients in need of a root canal treatment.
Clinical Benefit: Allows the dentist to locate and prepare the root canal for rotary instrumentation.

PRODUCT COMPOSITION

These instruments contain Stainless Steel, Polyetherimide, printing inks, and non-latex Silicone.

INDICATIONS FOR USE

Hand Files are instruments used for shaping, cleaning, and enlarging root canals for injured pulp conditions.

CONTRAINDICATIONS

None known.

WARNINGS

- Hand Files are non-sterile products and must be cleaned and sterilized prior to first use to prevent possible biohazard exposure.
- These products contain stainless steel and should not be used for individuals with known allergic sensitivity to this metal. The dental professional should inform the patient that stainless steel is known to cause allergic reactions.
- Use of a rubber dam is recommended when using these instruments for root canal procedures to avoid accidental aspiration or ingestion and avoid injury to the patient's oral soft tissues.

PRECAUTIONS

- Inspect the packaging and the instrument before use: if the packaging or the instrument are damaged, do not use and discard it.
- Determining the working length with the Endo Stop is imperative to ensure proper instrumentation.
- Ensure straight-line access and do not force a file down the canal as it might separate, get stuck, create transportation an/or ledge formation.
- Ensure that excessive stress is not placed on file during instrumentation as it might separate inside the root canal.
- Ensure that hand files are used in the correct sequence to avoid ledge formation, canal transportation or file separation.
- Ensure the final file used has a right size; it is not too large or too small to avoid excessive or insufficient instrumentation as it might cause clinical failure due to perforation, canal transportation or presence of infected pulp tissue left inside the root canal.
- Ensure the file is used in a root canal with proper irrigation/lubrification.
- Use personal protective equipment such as gloves and safety glasses.

ADVERSE REACTIONS

If a serious incident occurs with this medical device, report it to the legal manufacturer, Kerr Italia, S.r.l., and to the competent authority for that country in which the user and/or patient is established.

STEP-BY-STEP INSTRUCTIONS

Limitations on Use:

- Hand Files can be used a maximum of five (5) times.
- For extremely curved canals, use a new file for each canal.
- When the file shows: untwisted threads, dull cutting surfaces, bent instrument, signs of corrosion, loss of color marking and unable to sterilize the instrument prior to use.

Hand Files Technique

a. Obtain straight-line access.
b. Locate the canal orifices and obtain coronal patency.
c. Begin crown down by taking a coronal shaper to resistance.
d. Establish working length using hand files and apex locator. A radiograph may be taken to assist in working length determination.
e. In more difficult, curved, and narrow canals you must go through several recapitulations of the sequence.

Clinical Tips

- Irrigate the canal using NaOCl, 17% EDTA or a suitable endodontic lubricant agent to facilitate debris removal.
- Use light pressure and avoid force. When the instrument no longer advances apically, proceed to a smaller file.
- Prior to use, check instruments for defects; including deformation, breakage, fractures, corrosion, handle damage, loss of color coding or size marking. Discard any product with defects.
- If an instrument failes to function normally, it must not be used.

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